
Calf Implants



Calf Implants

Sientra calf implants are used in cosmetic surgery, as well as in reconstructive surgery to correct unevenness or irregularities of the calf, such as those caused by congenital malformation, neurological diseases, or trauma. They are composed of a silicone elastomer envelope and are available in two types: Symmetrical and Asymmetrical.




	DIME	REFERENCE		
A cm	B cm	C cm	VOL cc	
15.0	4.7	1.5	70 сс	412-150
19.7	5.4	1.9	125 сс	412-200
23.5	5.8	1.9	175 сс	412-230

Symmetrical Calf Implant

Asymmetrical Calf Implant

Smooth Surface

	DIMEN	REFERENCE		
A cm	B cm	C cm	VOL cc	
6.6	2.7	1.2	15 сс	402-065
9.9	3.1	1.5	30 сс	402-100
15.8	3.1	1.6	50 сс	402-160
19.1	3.2	1.8	70 сс	402-190
26.2	4.2	2.4	140 сс	402-260
26.9	3.6	1.9	105 сс	402-270



CALF IMPLANT

- 1) INCISION IN THE SKIN
- 2) INCISION OF THE APONEUROSIS OF THE MEDIAL GASTROCNEMIUS MUSCLES
- 3) IMPLANT
- 4) INNER POPLITEAL SCIATIC NERVE



Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.



PO Box 1490 | Santa Barbara, CA 93116-1490 | Tel: 805.562.3500 | Toll Free: 888.708.0808 Fax: 805.562.8401 | www.sientra.com

MDC-0020 R2